

THE SCIG ALTERNATIVE

AVAILABLE SCIG SETS

Code	Tubing	Lumen	Needle Length (mm)	Gauge	Priming Volume (ml)	Box Qty
SUB-104-G27	36"	Single	4	27	0.23	50
SUB-106	36"	Single	6	27	0.23	50
SUB-106-S	2"	Single	6	27	0.12	50
SUB-109-G24	36"	Single	9	24	0.23	50
SUB-109-G27	36"	Single	9	27	0.23	50
SUB-112-G24	36"	Single	12	24	0.23	50
SUB-204-G27	36"	Bi-furcated	4	27	0.40	10
SUB-212-G24	36"	Bi-furcated	12	24	0.40	10
SUB-212-G27	36"	Bi-furcated	12	27	0.40	10
SUB-250	36"	Bi-furcated	6	27	0.40	10
SUB-260	36"	Bi-furcated	9	27	0.40	10
SUB-270	3.5"	Bi-furcated	n/a	n/a	0.14	10
SUB-309-G24	36"	Tri-furcated	9	24	0.60	10
SUB-310	36"	Tri-furcated	6	27	0.65	10
SUB-312-G24	36"	Tri-furcated	12	24	0.60	10
SUB-312-G27	36"	Tri-furcated	12	27	0.60	10
SUB-320	36"	Tri-furcated	9	27	0.65	10
SUB-400	36"	Quad-furcated	6	27	0.80	10
SUB-409-G24	36"	Quad-furcated	9	24	0.80	10
SUB-410	36"	Quad-furcated	9	27	0.80	10
SUB-412-G24	36"	Quad-furcated	12	24	0.80	10
SUB-412-G27	36"	Quad-furcated	12	27	0.80	10
SUB-414-G27	36"	Quad-furcated	14	27	0.80	10
SUB-470	5"	Quad-furcated	n/a	n/a	0.30	10
SUB-506	36"	Penta-furcated	6	27	1.00	10
SUB-509	36"	Penta-furcated	9	27	1.00	10
SUB-512-G24	36"	Penta-furcated	12	24	1.00	10
SUB-606	36"	Hexa-furcated	6	27	1.20	10
SUB-609	36"	Hexa-furcated	9	27	1.20	10
SUB-612-G24	36"	Hexa-furcated	12	24	1.20	10

